



Information sheet PTS

ND virus detection – VLDIA328

These products are produced by the GD Animal Health and meant to be used in the Proficiency testing scheme.
(Spiked medium on FTA cards – Newcastle disease virus, NDV)

Samples

The sample set consists of 8 coded samples. The samples are spiked medium containing Newcastle disease virus. The samples have to be tested in duplicate (in different test runs) using the method(s) (of choice) of your laboratory. The samples should be tested under normal routine test conditions.

- 1 Set PTS samples.
1 set consists of 8 coded samples (spiked medium on FTA Cards) containing Newcastle disease virus.
- The samples should be stored at 2-8 °C and tested within 5 days.
- Sample will be sufficient to perform the required tests.

These products are produced, bottled and sealed at GD Animal Health and meant to be used in the Proficiency testing scheme.
No commercial value

All samples are treated with β -Propiolactone.

Reporting your findings

You can use Nucleic acid Amplification Techniques (NAT), e.g. PCR.
Results can be reported as positive-negative and/or Ct values.